

HD XRF Analyzer

for Screening & Quantification of Pb
in Toys



November 7, 2008



15 Tech Valley Drive
East Greenbush, N.Y. 12061
518.880.1500
www.xos.com





Company Highlights

Founded: 1990

Location: East Greenbush, NY

Products:

- Application specific X-Ray Analyzers
- X-Ray Optics & Beam Systems




better analysis counts

HDXRF



Demonstration

- Test #1: Quantification: 45 ppm Pb in 30 μ m paint on ABS plastic
- Test #2: Screening: 80 ppm Pb in 30 μ m paint on ABS plastic
- Test #3: Quantification: 80 ppm Pb in 30 μ m paint on ABS plastic
- Test #4: Quantification: 80 ppm Pb in 15 μ m paint on ABS plastic – 2 mm spot
- Test #5: Quantification: SRM 8115A 94 ppm in homogenous sample- 20s measurement
- Test #6: Screening: 160 ppm in 30 μ m paint on irregularly shaped wooden toy

Pb levels in the paint layers are correlated to ICP/OES and AA by outside test labs







Conventional XRF



- Unable to quantify or even detect low concentrations (poor signal-to-noise ratio)
- Average across multiple features and materials (large area illuminated)
- Average of the paint and the substrate

HD XRF




- Detect and quantify trace concentrations (excellent signal-to-noise) in:
 - paint
 - substrate
- Results for small features

Screening Mode

- Typical: 1 min. screening – user selectable
- 90 ppm level in 20um paint layer
- For thinner layers – extend screening time
- Low background enables fast screening for low concentrations of Pb in both paint and substrate
- Measurement averages results from substrate and paint layer
- Any Pb concentration > detection limit is reported as

Pb DETECTED

Low Noise → Increased Sensitivity




Example: Only 1 second of data collection

Test # 2: Screening: 80 ppm Pb in 30 um paint on ABS


Quantification Mode

- Quantifies and reports Pb concentration in:
 - the substrate
 - paint layers
- HD XRF: Multiple optic channels enable selective use of monochromatic excitation beams at selected energy levels
- Measurement time ~ 4 min.
- Low background ensures low detection limit for Pb in paint and substrate



Test #3: Quantification: 80 ppm Pb in 30 um paint on ABS

User Interface – Quantification Mode



Wide View

Spot View

Results

Element	Coating (ppm)	Base Material (ppm)
Pb	120	< 2
Sb	< 100	< 20
As	< 8	< 2
Ba	< 40	< 20
Cd	< 100	< 20
Cr	115	< 6
Hg	< 10	< 2
Se	< 8	< 2
Br	220	< 1
Cl	< 100	< 100

- Multi-element capacity
- User selects limits
- Encrypted secure data for recordkeeping and verification.
- Two images
 - Whole toy
 - Measurement spot

All images date & time stamped

Measuring Fine Details




HDX 1000
Analysis Area



Test #4: Quantification: 80 ppm Pb in 15 um paint on ABS – 2 mm spot


Small Spot Capability

- Small measurement area capability eliminates pass/fail ambiguity.
- $\mu\text{g}/\text{cm}^2$ vs. ppm (wt)
- Enables accurate and repeatable results of curved and irregular shaped surfaces.




Paint Layer Thickness Limitation

- Screening Mode:
 - In 1 min detect presence of Pb
 - Ensures concentration in the paint does not exceed 90 ppm in 20 μm paint layer
- Quantification Mode:
 - Quantify Pb at 80 ppm in 5 μm paint layer.
- Porous substrates:
 - substrate / paint layer interface not defined.



Precision and Accuracy

Limits of Detection

	Pb	Cd	Cr	As	Br	Sb	Se	Hg	Ba	Cl
Uniform PE	1	5*	3	0.8	0.6	2	0.8	1	3	50
Paint/Layers	8	40*	10	6	5	15	6	8	20	200

Unless noted all data is 3 sigma for 60 seconds

Extending measurements to 300 seconds would be half of these values

*Cadmium data reflects longer measurement time over 240 seconds.

$$L.L.D. = 3 \times k \times \sqrt{\frac{I_{back}}{T}}$$

Test #5: Quantification: SRM 8115A 94 ppm in homogenous sample



30 μ m paint on
ABS plastic


Repeatability

ABS-L-T 160 ppm nominal	ABS-UL-T 80 ppm nominal	NMIJ8115 94 ppm
166	94	97
150	90	96
157	89	95
157	84	97
143	99	98
166	87	96
183	77	95
162	80	98
155	86	98
169	74	97
149	78	95
166	88	97
141	73	94
157	87	96
154	80	98
160	78	96
164	78	95
149	93	96
162	80	96
168	93	98


	ABS-L-T	ABS-UL-T	NMIJ8115
min	141	73	94
average	159	84	96.4
mode	166	78	96
median	159	85	96
max	183	99	98
StDev	10	7	1.23
StDev %	6%	9%	1%

Accuracy

Paint Layers and Substrates



30 μ m paint on
ABS plastic



Comparative Data


Lead Levels in 30 micron Paint on ABS Plastic Substrate				
HD XRF	60	143	790	1301
FAAS	66	126	750	1296
ICP-OES	79	151	847	1400

- All data in ppm
- Samples prepared by XOS
- FAA conducted by Wadsworth Laboratories
- ICP-OES conducted by Galbraith Laboratories


Irregular Shapes



Conventional XRF



HD XRF



Test #6: Screening: 160 ppm in 30 um paint on curved or irregular surface

Substrate Matrices

- Plastics
- Wood
- Metal
- Leather
- Fabric
- Paper / Cardboard
- Glass and Ceramics



Summary of Benefits & Features

- Quantification of both Paint Layers & Substrates
 - LOD for Pb:
 - Paint 8 ppm
 - Homogeneous sample 1 ppm (PE substrate)
- Screening at 90 ppm in paint and substrate
- 10 Elements (incl. ASTM 963, + Br and Cl)
- Small Analysis Area: 1 mm
- Accurate and repeatable results with irregular shapes
- Variety of Substrate Materials

HDX 1000 Analyzer



in Toys Analyzer 玩具铅含量分析仪



Factories ▪ Test Labs ▪ Regulatory ▪ Retail

⁸² Pb Lead	⁴⁸ Cd Cadmium	⁵¹ Sb Antimony	³⁵ Br Bromine	³³ As Arsenic	⁸⁰ Hg Mercury	²⁴ Cr Chromium	³⁴ Se Selenium	⁵⁶ Ba Barium	¹⁷ Cl Chlorine
------------------------------------	---------------------------------------	--	---------------------------------------	---------------------------------------	---------------------------------------	--	--	--------------------------------------	--

工厂 • 实验室 • 出口检验 • 零售商