



AN ENVIRONMENTAL ANALYTICAL LABORATORY

COMPREHENSIVE VALIDATION PACKAGE

ATL Applications

INVENTORY SHEET

WORK ORDER # 1010262A

	Page Nos.	
	From	To
1. Work Order Cover Page & Laboratory Narrative & Table	1	2
2. Sample Results and Raw Data (Organized By Sample)	3	6
a. ATL Sample Results Form		
b. Target Compound Raw Data		
-Internal Standard Area and Retention Time Summary (If Applicable)		
-Surrogate Recovery Summary (If Applicable)		
-Chromatogram(s) and Ion Profiles (If Applicable)		
3. QC Results and Raw Data		
a. Method Blank (Results + Raw Data)	-	-
b. Surrogate Recovery Summary Form (If Applicable)	-	-
c. Internal Standard Summary Form (If Applicable)	-	-
d. Duplicate Results Summary Sheet	-	-
e. Matrix Spike/Matrix Spike Duplicate (Results + Raw Data)	-	-
f. Initial Calibration Data (Summary Sheet + Raw Data)	-	-
g. MDL Study (If Applicable)	-	-
h. Continuing Calibration Verification Data	-	-
i. Second Source LCS (Summary + Raw Data)	-	-
j. Extraction Logs	-	-
k. Instrument Run Logs/Software Verification	7	14
l. GC/MS Tune (Results + Raw Data)	-	-
4. Shipping/Receiving Documents:		
a. Login Receipt Summary Sheet	15	16
b. Chain-of-Custody Records	17	17
c. Sample Log-In Sheet	18	18
d. Misc. Shipping/Receiving Records (list individual records)		
<u>Sample Receipt Discrepancy Report</u>	19	20
5. Other Records (describe or list)		
a. <u>Manual Spectral Defense</u>	-	-
b. <u>Manual Intergrations</u>	-	-
c. <u>Manual Calculations</u>	-	-
d. <u>Canister Dilution Factors</u>	-	-
e. <u>Laboratory Corrective Action Request</u>	-	-
f. <u>CAS Number Reference</u>	21	22
g. <u>Variance Table</u>	-	-
h. <u>Canister Certification</u>	-	-
i. <u>Data Review Check Sheet</u>	23	23

Completed by:

V. Belitsky

(Signature)

Vera Belitsky/ Document Control

(Print Name & Title)

10/27/10

(Date)


WORK ORDER #: 1010262A

Work Order Summary

CLIENT:	Mr. Brian Baker Environmental Health & Engineering, Inc. 117 Fourth Avenue Needham, MA 02494	BILL TO:	Accounts Payable Environmental Health & Engineering, Inc. 117 Fourth Avenue Needham, MA 02494
PHONE:	800-825-5343	P.O. #	17314
FAX:	781-247-4305	PROJECT #	17314
DATE RECEIVED:	10/13/2010	CONTACT:	Ausha Scott
DATE COMPLETED:	10/26/2010		

<u>FRACTION #</u>	<u>NAME</u>	<u>TEST</u>
01A	116171	ATL Applications
02A	116172	ATL Applications
03A	116173	ATL Applications
04A	116174	ATL Applications
05A	116175	ATL Applications
06A	116176	ATL Applications
07A	Lab Blank	ATL Applications
07B	Lab Blank	ATL Applications
08A	LCS	ATL Applications

CERTIFIED BY:



Laboratory Director

DATE: 10/26/10

This report shall not be reproduced, except in full, without the written approval of Air Toxics Ltd.

180 BLUE RAVINE ROAD, SUITE B FOLSOM, CA - 95630
(916) 985-1000 . (800) 985-5955 . FAX (916) 985-1020

**LABORATORY NARRATIVE
Hydrogen Sulfide by Radiello 170
Environmental Health & Engineering, Inc.
Workorder# 1010262A**

Six Radiello 170 (H₂S) samples were received on October 13, 2010. The procedure involves adsorption of H₂S by zinc acetate to form zinc sulfide. The sulfide is then recovered by extraction with water and addition of ferric chloride in a strongly acidic solution to produce methylene blue. Methylene blue absorbance is then measured at 665 nm using a spectrophotometer. Results are reported in uG and uG/m³.

Sampling rate of 69 mL/min for H₂S was provided by the manufacturer.

Receiving Notes

Sample collection date was not provided on the Chain of Custody for samples 116171, 116172, 116173, 116174, 116175 and 116176. The client was contacted and a date of 9/28/10 was provided.

Analytical Notes

Results were calculated based on 25 deg C without temperature correction. The actual exposure time was used to calculate sample concentrations and reporting limits.

An exposure time of 19805 minutes was used for the QC samples and trip blanks.

Definition of Data Qualifying Flags

Eight qualifiers may have been used on the data analysis sheets and indicate as follows:

- B - Compound present in laboratory blank greater than reporting limit.
- J - Estimated value.
- E - Exceeds instrument calibration range.
- S - Saturated peak.
- Q - Exceeds quality control limits.
- U - Compound analyzed for but not detected above the detection limit.
- M - Reported value may be biased due to apparent matrix interferences.
- N - The identification is based on presumptive evidence.

File extensions may have been used on the data analysis sheets and indicates as follows:

- a-File was requantified
- b-File was quantified by a second column and detector
- r1-File was requantified for the purpose of reissue

Sample Results and Raw Data

AIR TOXICS LTD.

ATL Application # 59 for RAD 170 (Hydrogen Sulfide)

Spectrophotometer

Field Sample I.D.	Lab Sample I.D.	Collection Date	Analysis Date	Dilution Factor	Reporting Limit (ug)	Reporting Limit (ug/m3)	Amount (ug)	Amount (ug/m3)
116171	1010262A-01A	NA	10/18/2010	1.00	0.80	0.55	ND	ND
116172	1010262A-02A	NA	10/18/2010	1.00	0.80	0.55	ND	ND
116173	1010262A-03A	NA	10/18/2010	1.00	0.80	0.55	ND	ND
116174	1010262A-04A	NA	10/18/2010	1.00	0.80	0.55	ND	ND
116175	1010262A-05A	NA	10/18/2010	1.00	0.80	0.55	ND	ND
116176	1010262A-06A	NA	10/18/2010	1.00	0.80	0.55	ND	ND
Method Blank	1010262A-07A	NA	10/18/2010	1.00	0.80	0.55	ND	ND
Method Blank	1010262A-07B	NA	10/18/2010	1.00	0.80	0.55	ND	ND
LCS	1010262A-08A	NA	10/18/2010	1.00	0.80	0.55	<div>%Rec</div> <div>128</div>	

COMMENTS: 1. NA=Not Applicable

2. ND=Not Detected

3. Exposure time of 19805 minutes was assumed for the QC samples.

4. Background subtraction not performed.

Hydrogen Sulfide Radiello Calculation Worksheet

Workorder #: 1010262A

Sampling Rate (ng/ppb.min)

0.096 Typically 0.096 for H2S

Sampling T (deg C)

2.5 Typically 25

Volume (mL)

10.5 Typically 10.5 for H2S

Date of Analysis:

8/2010

Corrected Q

0.096

Takes into account temp

LabSampleID

Client

Date of Collection **A**

Abs

Duration (min)

DF

Conc (ug/mL) of sulfide

Conc (ug) of sulfide

Conc (ug) of H2S

Conc (ppb) of H₂SConc ($\mu\text{g}/\text{m}^3$) of H_2S

Q includes conversion from Sulfide to H₂S

$$\text{Conc (}\mu\text{g)} \times 1000$$

Q x Duration

ppbx mw

24.45

T Corrected, no Blank correction

[illegible]

QC Duration

Duration

CCV Spike Amt

the Art

Verified: HH and AW on 9/4/05

QC Results and Raw Data

Work Order: 1010262A/100269D

Date: 10/18/10

Method: Rad 170

Analyst: M. SKIsmore


Wavelength: 665nm

Standard ID	Concentration	ABS
	sulfide (mg/mL)	
Level 1 1993-80-E	0.0716	0.097
Level 2 -D	0.143	0.180
Level 3 -C	0.286	0.356
Level 4 -B	0.572	0.683
Level 5 -A	1.145	1.237
ICV 1993-81	0.286	0.345

$$r = \frac{0.9974}{1.062}$$

$$m = \frac{0.0396}{0.0396}$$


ICV % Recovery = 101

Fraction	Dilution	ABS ^{MJS 10/18/10}	Sample ID	Sample Volume	Comments
01A	1.00	0.097 0.040	116171	10.5 mL	
02A		0.118 0.039	116172		
03A		0.050	116173		
04A		0.035	116174		
05A		0.021	116175		
06A		0.021	116176		
49A		0.091	115813		
50A		0.118	115814		
51A		0.123	^{refilled MJS} 115815		
52A		0.131	115816		
53A		0.029	115817		
54A		0.026	115818		
54AA		0.025			
B1K1		0.026	N/A		Lot: 10101
B1K2		0.027			
LC5		0.221			0.133 mg/L
CCV		0.361			0.286 mg/mL
					
MJS 10/19/10					

Procedure:

- 1.) Add 10 mL of H₂O to sample tube, cap and vortex for 1 minute.
- 2.) Add 0.5 mL of Ferric Chloride-Amine solution and cap immediately.
- 3.) Allow color to develop for 30 minutes.
- 4.) Measure absorbance at 665nm.

MJS 10/19/10


 Signed

 10/19/10
 Date

Spectrophotometer Standard Preparation Log

@Air Toxics Ltd. Log Book #: 1993

Standard ID: 1993-80

Project: Rad 170 calibration Curve

Analyst: M. Skidmore

Preparation Date: 10/18/10

Expiration Date: 10/18/10

Solvent: HPLC H₂O

Solvent Lot #: DB 270

Procedure/Comments:

Solution A: 2 mL of Code Rad 171 (1476-1736, exp 2/3/11) (located in ER1B) with 98 mL of D.I. H₂O = 1.145 µg/mL

Solution B: 2.5 mL of Solution A with 2.5 mL of D.I. H₂O = 0.572 µg/mL


Solution C: 1.25 mL of Solution A with 3.75 mL of D.I. H₂O = 0.286 µg/mL

Solution D: 0.625 mL of Solution A with 4.375 mL of D.I. H₂O = 0.143 µg/mL

Solution E: 0.375 mL of Solution A with 5.625 mL of D.I. H₂O = 0.0716 µg/mL

Note: Each solution was measured immediately after it was prepared. Solution A is only stable in the flask it was prepared in.

MJS 10/18/10



MJS 10/18/10

Spectrophotometer Standard Preparation Log

@Air Toxics Ltd. Log Book #: 1993

Standard ID: 1993-81 *MTS 10/18/10*
Project: Rad 170 *ICV*
Analyst: *FM*
Preparation Date: 10/18/10
Expiration Date: 10/18/10

Solvent: *HPLC water*
Solvent Lot #: *DB270*


Procedure/Comments: _____

Solution A: 2 mL of Code Rad 171 (1476-1736, exp 2/3/11) (located in ER1B) with
98 mL of D.I. H₂O = 1.145 µg/mL

Solution C: 1.25 mL of Solution A with 3.75 mL of D.I. H₂O = 0.286 µg/mL

Note: Each solution was measured immediately after it was prepared. Solution A is only
stable in the flask it was prepared in.

MTS 10/18/10



MTS 10/18/10

Fauzin
Signed

10/18/10
Date

[Signature]
Reviewed

10/18/10
Date

Spectrophotometer Standard Preparation Log

@Air Toxics Ltd. Log Book #: 1993Standard ID: 1993-79Project: Rad 170 H₂S LCSAnalyst: M. SkidmorePreparation Date: 10/18/10Expiration Date: 10/18/10Solvent: HPLC H₂OSolvent Lot #: DB 270

Procedure/Comments: _____

A Rad 170 cartridge (lot: 10101) was placed in a 40 mL VOA vial. 10.0 mL of D.I. H₂O was aliquoted into the vial. 1.0 mL of H₂S gas (1476-1497; 1000 ppm) was injected into the vial, into the H₂O. The solution was allowed to gently shake for 2 hours. Then 0.5 of the ferric-chloride-amine (1993-78) was added to the vial and capped immediately. The solution was allowed to sit for 30 minutes and the absorbance was measured at 665 nm.

MJS 10/18/10

MJS
10/18/10

Spectrophotometer Standard Preparation Log

@Air Toxics Ltd. Log Book #: 1993

Standard ID: 1993-78

Project: Ferric Chloride-Amine Solution Rad170

Analyst: M. Skidmore


Preparation Date: 10/18/10

Expiration Date: 10/18/10

Solvent: HPLC H₂O

Solvent Lot #: DB270

Procedure/Comments: Add 12.5 mL of ferric chloride solution
(1993-77, exp 10/18/11) with 62.5 mL of amine solution
(1993-76, exp 11/18/10).



MS
10/18/10

Spectrophotometer Standard Preparation Log

@Air Toxics Ltd. Log Book #: 1993

Standard ID: 1993-76
Project: Rad 170 Amine Solution
Analyst: M. Skidmore
Preparation Date: 10/18/10
Expiration Date: 11/18/10

Solvent: HPLC H₂O
Solvent Lot #: DB 270

Procedure/Comments:

Sulfuric Acid Solution:

Slowly add 6.25 mL of concentrated sulfuric acid to 2.5 mL of D.I. H₂O, and let the solution cool. (sulfuric acid lot: 01428LS).

Amine Solution:

Dissolve 1.6875g of N,N-dimethyl-p-phenyldiammonium oxalate (located in ER1A; Lot: 63797PJ) in the above mentioned sulfuric acid solution. Dilute this solution to 250 mL with sulfuric acid-water 1:1 v/v. (This is roughly 120 mL H₂O + 120 mL sulfuric acid).

MJS 10/18/10

MJS 10/18/10

Spectrophotometer Standard Preparation Log

@Air Toxics Ltd. Log Book #: 1993

Standard ID: 1993-77

Project: Ferric Chloride Solution Rad 170

Analyst: M. Skidmore


Preparation Date: 10/18/10

Expiration Date: 10/18/11

Solvent: HPLC H₂O

Solvent Lot #: DB 270

Procedure/Comments: Dissolve 125 g of ferric chloride hexahydrate
(located in ERAC, lot: 73297) in 50 mL of H₂O,



MJS 10/18/10

Shipping/ Receiving Documents

180 Blue Ravine Road, Suite B
Folsom, CA 95630

Phone (916) 985-1000 FAX (916) 985-1020
Hours 8:00 A.M. to 6:00 P.M. Pacific

COMPANY: Environmental Health & Engineering, Inc.
ATTENTION: Mr. Brian Baker
FAX #: 781-247-4305
FROM: Sample Receiving
Workorder #: 1010262A
of pages (Including Cover): 4
10/27/2010

Thank you for selecting Air Toxics Ltd. We have received your samples and have found discrepancies. In order to expedite analysis and reporting, please review the attached information for accuracy.

Corrections can be faxed to **Ausha Scott at 916-985-1020.**

ATL will proceed with the analysis as specified on the Chain of Custody and Sample Login page.

In accordance with your company's contract, this account is required to have a PO that is fully executed by both parties which also covers the cost of the workorder before any data can be released. Please ensure that you have given all appropriate information to our Project Manager so that there will be no delay in reporting of the data you are requesting.

The following discrepancy has been observed:

Samples were received without documentation regarding collection date on the Chain of Custody. The sampling date of 9/28/10 you have provided by telephone/fax/e-mail will be used to determine the extent of hold time.

Your prompt response is appreciated.

1010262

Environmental
Health &
Engineering, Inc.

CHAIN OF CUSTODY FORM

DATE: 12 OCT 10

FROM: Environmental Health and Engineering, Inc.
117 Fourth Avenue
Needham, MA 02494-2725

TO: AIR TOXICS

Please send invoices to ATTN: Accounts Payable
Please send reports to ATTN: Data Coordinator

In all correspondence regarding this matter, please refer to EH&E Project # 17314The cost of this analysis will be covered by EH&E Purchase Order # 17314For EH & E Data Coordinator - URGENT DATA ☒

SAMPLE ID	SAMPLE TYPE	ANALYTICAL METHOD/NUMBER	OTHER: Time/Date/Vol.
01A 116171	AIR / PASSIVE	H ₂ S ANALYSIS	13D 18H 5M
02A 116172	I	I	I
03A 116173	I	I	I
04A 116174	I	I	I
05A 116175	I	I	ØD
06A 116176	I	I	I
116177	AIR / PASSIVE	FORMALDEHYDE ANALYSIS	13D 18H 5M
116178	I	I	I
116179	I	I	I
116180	I	I	I
116181	I	I	ØD

Special instructions:

☒ Standard turn around time☐ Rush by _____ date/time☐ Other _____☐ Fax results 781-247-4305☐ RETURN SAMPLES☒ Electronic transfer - datacoordinator@ehinc.com☒ Additional report recipient blakere@ehinc.com

Each signatory please return one copy of this form to the above address

Relinquished by: [Signature] of Environmental Health & Engineering, Inc.Date: 10/12/10Received by: [Signature] of (company name) ATLDate: 10/13/10 09:00

Relinquished by: _____ of (company name) _____

Date: _____

Received by: _____ of (company name) _____

Date: _____

Relinquished by: _____ of (company name) _____

Date: _____

Received by: _____ of (company name) _____

Date: _____

Lab Data

Received by: _____ of Environmental Health & Engineering, Inc.

Date: _____



FedEx

Page 1 of 1

SAMPLE RECEIPT SUMMARY

WORKORDER 1010262A

Client	Phone	Date Promised: 10/26/10 11:59 pm
Mr. Brian Baker	800-825-5343	Date Completed:
Environmental Health & Engineering, Inc.	Fax	Date Received: 10/13/10
117 Fourth Avenue	781-247-4305	PO#: 17314
Needham, MA 02494		Project#: 17314
Sales Rep: TL		Total \$: \$ 510.00
		Logged By: AW

<u>Fraction</u>	<u>Sample #</u>	<u>Analysis</u>	<u>Collected</u>	<u>Amount\$</u>
01A	116171	ATL Applications	NA	\$80.00
02A	116172	ATL Applications	NA	\$80.00
03A	116173	ATL Applications	NA	\$80.00
04A	116174	ATL Applications	NA	\$80.00
05A	116175	ATL Applications	NA	\$80.00
06A	116176	ATL Applications	NA	\$80.00
Misc. Charges eCVP (6) @ \$5.00 each.				\$30.00

Note: Samples received after 3 P.M. PST are considered to be received on the following work day.
Atlas Project Name/Profile#: CPSC/14482

BILL TO: Accounts Payable
Environmental Health & Engineering, Inc.
117 Fourth Avenue
Needham, MA 02494

Analysis Code: Other GC

TERMS:

Reporting Method: ATL Application #59 H2S-Radiello 170

180 BLUE RAVINE ROAD, SUITE B FOLSOM, CA - 95630
(916) 985-1000 . (800) 985-5955 . FAX (916) 985-1020

@ Air Toxics Ltd	Title: Sample Discrepancy Report		Release Date: 03/03/10
	Form #: F1.3	Revision #: 1	Revision Date: 10/7/08
			Page #: 1 of 2

Sample Discrepancy Report

Identification

Initiated By: AW Project ID: 14482 PM: AS Date: 10/13/2010 Discrepancy Type: ☐ 1. ☒ 2. ☐ 3.

Workorder(s) affected: 1010262A/B Sample(s) affected: All

1. Sample Receipt Discrepancies

Narration Not Required:

- 1.1. ☐ Sample container (cartridge/tube/VOA vial) was received broken, however sample was intact.
- 1.2. ☐ No brass cap on canister.
- 1.3. ☐ Date of Collection noted on first sample, but no arrow down to indicate all samples.

Notify Lab for further determination:

- 1.4. ☐ Tedlar bag received with minimal volume.

Initials: _____ Date: _____

Narration Required in Lab Narrative and Sample Confirmation:

- 1.5. ☐ COC was not filled out in ink.
- 1.6. ☐ COC improperly relinquished / received.
- 1.7. ☐ Sample tags / can numbers do not match the COC.
- 1.8. ☐ Sample date ☐ error / ☐ missing on COC but noted on sample tag (check one).
- 1.9. ☐ Custody Seal on the outside of the container was ☐ broken / ☐ improperly placed (check one).
- 1.10. ☐ ID-none on the sample Tag/Blank
- 1.11. ☐ Other (describe below).

Describe the Discrepancy: _____

2. Sample Receipt/Screening Discrepancies requiring PM notification

Document on Cover Page of Sample Receipt Confirmation and in Receiving Notes of Lab Narrative

If Section II. is filled out PM must be notified within 24 hrs of initiation

- | | |
|---|---|
| <ul style="list-style-type: none"> 2.1. <input type="checkbox"/> COC was not received with samples. 2.2. <input type="checkbox"/> Analysis method(s) is <input type="checkbox"/> not specified / <input type="checkbox"/> incorrectly specified (check one) on the COC. 2.3. <input type="checkbox"/> Incorrect sampling media / container for analysis requested. 2.4. <input type="checkbox"/> Number of samples on the COC does not match the number of samples that were received. 2.5. <input type="checkbox"/> Samples were received expired. 2.6. <input checked="" type="checkbox"/> Sampling date (time for sulfur) is not documented for <input type="checkbox"/> <u>some</u> / <input checked="" type="checkbox"/> <u>any</u> samples (check one). 2.7. <input type="checkbox"/> Sample received with amount of H₂O in the Tedlar Bag. 2.8. <input type="checkbox"/> Sample cannot be analyzed. Container was <input type="checkbox"/> received broken / <input type="checkbox"/> leaking / <input type="checkbox"/> flat / <input type="checkbox"/> defective. 2.9. <input type="checkbox"/> Tedlar bag / canister received emitting a strong odor; Sample <input type="checkbox"/> can / <input type="checkbox"/> cannot (check one) be analyzed. 2.10. <input type="checkbox"/> Tedlar Bag for Sulfur analysis has metal fitting. 2.11. <input type="checkbox"/> Environmental Supply Company valves 2.12. <input type="checkbox"/> Sorbent samples-sampling volume was not provided | <ul style="list-style-type: none"> 2.13. <input type="checkbox"/> Flow controller used – canister samples received at ambient or under pressure. 2.14. <input type="checkbox"/> Canister was at ambient pressure at time of pressurization and (check all that apply):
 <input type="checkbox"/> Canister failed leak check on two manifolds,
 <input type="checkbox"/> Canister valve was open,
 <input type="checkbox"/> Brass nut was loose/not present.
 <input type="checkbox"/> Sample can be analyzed
 <input type="checkbox"/> Cannot be analyzed 2.15. <input type="checkbox"/> Canister sample received with a vacuum difference >5.0"Hg between the receipt vac. And the final vac. reported on the COC, indicating loss of vacuum. 2.16. <input type="checkbox"/> Canister sample received at >15"Hg (<u>not</u> identified as a Trip/Field Blank). 2.17. <input type="checkbox"/> Canister Trip Blank received at low vacuum (< 25"Hg). 2.18. <input type="checkbox"/> Sorbent Sample received outside method required temperature of 2°C to 6°C; <input type="checkbox"/> ice / <input type="checkbox"/> blue ice (check one) was present. A temp. Blank <input type="checkbox"/> was / <input type="checkbox"/> was not present (check one). 2.19. <input type="checkbox"/> Other (describe below) |
|---|---|

Initials: _____ Date: _____ Notify Receiving: ☐ Notify PM: ☐

Describe the Discrepancy: _____

3. Lab Discrepancies requiring Team Leader/PM notification

Document in Analytical Notes of Lab Narrative

If Section III. is filled out PM must be notified within 24 hrs of initiation

- | | |
|--|--|
| 3.1. <input type="checkbox"/> Tedlar Bag found to be leaking at the time of analysis; sample <input type="checkbox"/> can / <input type="checkbox"/> cannot (check one) be analyzed. | 3.6. <input type="checkbox"/> Sample loss due to instrument malfunction / broken glassware. |
| 3.2. <input type="checkbox"/> Tedlar Bag found to be flat/low volume; sample cannot be analyzed. | 3.7. <input type="checkbox"/> Low/high surrogate recoveries noted in QC/sample(s) for extractable samples. |
| 3.3. <input type="checkbox"/> Sulfur samples received with insufficient time to analyze prior to expiration. | 3.8. <input type="checkbox"/> Reporting Limit was raised. |
| 3.4. <input type="checkbox"/> Canister found to be leaking at the time of analysis. | 3.9. <input type="checkbox"/> Post weight > Pre weight in field/lab Blank for PM10/TSP samples. |
| 3.5. <input type="checkbox"/> VOST tube saturated; bag dilution necessary. | 3.10. <input type="checkbox"/> Other (describe below). |

Initials: _____ Date: _____ Notify Receiving: ☐ Notify PM: ☐

Team Lead Initials: _____ Date: _____

Describe the Discrepancy: _____

How Does this Affect Client: _____

Project Manager Use Only

Project Manager Notification

☒ Section 2 Complete

☐ Section 3 Complete

Action:

- ☐ It is not necessary to notify the client. Narrate the discrepancy in Receiving Notes/Analytical Notes of Lab Narrative.

PM Initials: _____ Date: _____

- ☒ Client notification required. See attached client contact / email, or comments below:

Client Notification:

PM Initials: AS Person notified: B.Baker

Date: 10/13/2010

- ☐ Waiting for Client Reply

Comments: Client emailed spreadsheet on 10/18

☐ Notify Lab Name: _____ Date: _____ Notify Receiving: ☐

- ☐ Additional notifications attached.

Additional Comments:

Other Records



Method : ATL Application #59 H2S-Radiello 170

CAS Number	Compound	Rpt. Limit (ug)
7783-06-4	Hydrogen Sulfide	1.2

@ Air Toxics Ltd	Title: Data Review Checklist		Release Date: 07/28/10	
	Form #: F1.27	Revision #: 2	Revision Date: 07/27/10	Page #: 1 of 2

DATA REVIEW CHECKLIST

Work Order #:

1010262A

A₁ A₂ W T R Q

- ☒ ☐ ☒ ☐ ☐ ☐ Analysis/Reporting vs. Project Profile/SOP requirements checked (i.e. 100% Dups, J-Flag to MDL, etc)
- ☒ ☐ ☒ ☐ ☐ ☐ The final report has the correct reporting list, special units, and header info.
- ☒ ☐ ☒ ☐ ☐ ☐ Non-Standard sublist printed/verified, LOQ and LOD verified
- ☒ ☐ ☒ ☐ ☐ ☐ Lab Narrative is correct (proper method & description/Receiving & Analytical notes correct)
- ☒ ☐ ☒ ☐ ☐ ☐ Sample Discrepancy Report (SDR) is completed

- ☒ ☐ ☐ ☐ ☐ ☐ Corrective Action issued - # _____
- ☒ ☐ ☐ ☐ ☐ ☐ Unusual circumstances have been documented in the notes section below

LUMEN validation report present and initialed

CIRCLE (YES / NO)

- ☒ ☐ ☒ ☐ ☐ ☐ Lab Blank, CCV, LCS and DUP met QC criteria
- ☒ ☐ ☒ ☐ ☐ ☐ Hold time is met for all samples
- ☒ ☐ ☒ ☐ ☐ ☐ Appropriate data qualifier flags are applied
- ☒ ☐ ☒ ☐ ☐ ☐ Manual integrations for samples and QC are properly documented
- ☒ ☐ ☒ ☐ ☐ ☐ Samples analyzed within the project or method specific clock
- ☒ ☐ ☒ ☐ ☐ ☐ Retention times have been verified
- ☒ ☐ ☒ ☐ ☐ ☐ Appropriate ICAL(s) included, %RSD Recalculation

- ☒ ☐ ☒ ☐ ☐ ☐ At least one result per sample is verified against the target quant sheets/raw data
- ☒ ☐ ☒ ☐ ☐ ☐ Dilution factor correctly calculated (sample load volume, syringe and bag dilutions, can pressurization(s))
- ☒ ☐ ☒ ☐ ☐ ☐ Correct amount of sample analyzed (i.e. sample not over-diluted)
- ☒ ☐ ☒ ☐ ☐ ☐ Spectra verified - documentation of spectral defense included (Section 5A of eCVP pkg)

- ☐ ☐ ☐ ☐ ☐ ☐ TICs resemble reference spectra
- ☐ ☐ ☐ ☐ ☐ ☐ TICs between duplicate samples are consistent
- ☒ ☐ ☒ ☐ ☐ ☐ Checked samples for trends (i.e. Influent vs. Effluent, Field Dups, Field/Trip Blank, etc.)
- ☒ ☐ ☒ ☐ ☐ ☐ Data for multiple analyses of sample(s) has been evaluated for comparability of results

- ☐ ☐ ☒ ☐ ☐ ☐ Special units for all samples in the final report are correctly calculated
- ☐ ☐ ☒ ☐ ☐ ☐ Manually entered results checked (i.e. TPH/NMOC)

- ☒ ☐ ☒ ☐ ☐ ☐ Chain of Custody verified for any special comments (i.e. different compounds/RLs, action levels)
- ☐ ☐ ☒ ☐ ☐ ☐ Chain of Custody scanned correctly
- ☒ ☐ ☒ ☐ ☐ ☐ Verify sample id's vs. chain of custody
- ☒ ☐ ☒ ☐ ☐ ☐ Date MDL(s) performed per instrument(s)

- ☒ ☐ ☒ ☐ ☐ ☐ Samples pressurized w/ appropriate gas (N₂ or He) ☐ Other (i.e. Tedlar bag, cartridge, sorbent)
- ☐ ☐ ☐ ☐ ☐ ☐ Final pressure consistent with canister size (6L vs. 1L)
- ☐ ☐ ☐ ☐ ☐ ☐ Verify receipt pressures

- ☒ ☐ ☐ ☐ ☐ ☐ Verify canister ID #'s
- ☐ ☐ ☐ ☐ ☐ ☐ Final invoice amount correct (adjusted for TAT, Penalties, Re-issue Charges etc.)
- ☒ ☐ ☐ ☐ ☐ ☐ Final PDF report reviewed for correctness

Notes: (to include: noting samples with QA/QC problems, Blanks with positive hits, narratives, etc.)

A/R: 19,805 minutes duration used for all QC's and Trip Blanks,

T/Q:

A ₁ /A ₂ (Analytical Review/Date)	W/T (Write-up/Tech Review/Date)	R* (Report Review/Date)	Q (QA Review/Date)
A ₁ : <u>Milano 10/27/10</u>	W: <u>Milano 10/27/10</u>	R: _____	_____

A₂:

T:

Note (1): Please check all the appropriate boxes. Indicate "NA" for any statement that does not apply.

Note (2): Report reviewer and write-up reviewer must be separate individuals for DoD & Client Specific projects.

* Report Review is completed for DoD & Client Specific projects only.