



AN ENVIRONMENTAL ANALYTICAL LABORATORY

COMPREHENSIVE VALIDATION PACKAGE

ATL Applications

INVENTORY SHEET

WORK ORDER # 1011113A

	Page Nos.	
	From	To
1. Work Order Cover Page & Laboratory Narrative & Table	1	3
2. Sample Results and Raw Data (Organized By Sample)	4	7
a. ATL Sample Results Form		
b. Target Compound Raw Data		
-Internal Standard Area and Retention Time Summary (If Applicable)		
-Surrogate Recovery Summary (If Applicable)		
-Chromatogram(s) and Ion Profiles (If Applicable)		
3. QC Results and Raw Data		
a. Method Blank (Results + Raw Data)	-	-
b. Surrogate Recovery Summary Form (If Applicable)	-	-
c. Internal Standard Summary Form (If Applicable)	-	-
d. Duplicate Results Summary Sheet	-	-
e. Matrix Spike/Matrix Spike Duplicate (Results + Raw Data)	-	-
f. Initial Calibration Data (Summary Sheet + Raw Data)	-	-
g. MDL Study (If Applicable)	-	-
h. Continuing Calibration Verification Data	-	-
i. Second Source LCS (Summary + Raw Data)	-	-
j. Extraction Logs	-	-
k. Instrument Run Logs/Software Verification	8	15
l. GC/MS Tune (Results + Raw Data)	-	-
4. Shipping/Receiving Documents:		
a. Login Receipt Summary Sheet	16	17
b. Chain-of-Custody Records	18	18
c. Sample Log-In Sheet	19	20
d. Misc. Shipping/Receiving Records (list individual records)		
<u>Sample Receipt Discrepancy Report</u>	21	22
5. Other Records (describe or list)		
a. <u>Manual Spectral Defense</u>	-	-
b. <u>Manual Intergrations</u>	-	-
c. <u>Manual Calculations</u>	-	-
d. <u>Canister Dilution Factors</u>	-	-
e. <u>Laboratory Corrective Action Request</u>	-	-
f. <u>CAS Number Reference</u>	23	24
g. <u>Variance Table</u>	-	-
h. <u>Canister Certification</u>	-	-
i. <u>Data Review Check Sheet</u>	25	25

Completed by:

Kara McKiernan

(Signature)

Kara McKiernan/ Document Control

(Print Name & Title)

11/16/10

(Date)

WORK ORDER #: 1011113A

Work Order Summary

CLIENT:	Mr. Brian Baker Environmental Health & Engineering, Inc. 117 Fourth Avenue Needham, MA 02494	BILL TO:	Accounts Payable Environmental Health & Engineering, Inc. 117 Fourth Avenue Needham, MA 02494
PHONE:	800-825-5343	P.O. #	17131
FAX:	781-247-4305	PROJECT #	17131
DATE RECEIVED:	11/04/2010	CONTACT:	Ausha Scott
DATE COMPLETED:	11/15/2010		

<u>FRACTION #</u>	<u>NAME</u>	<u>TEST</u>
01A	118581	ATL Applications
02A	118582	ATL Applications
03A	118583	ATL Applications
04A	118584	ATL Applications
05A	118585	ATL Applications
06A	118586	ATL Applications
07A	118597	ATL Applications
08A	118598	ATL Applications
09A	118599	ATL Applications
10A	118600	ATL Applications
11A	118601	ATL Applications
12A	118602	ATL Applications
13A	118613	ATL Applications
14A	118614	ATL Applications
15A	118615	ATL Applications
16A	118616	ATL Applications
16AA	118616 Lab Duplicate	ATL Applications
17A	Lab Blank	ATL Applications

Continued on next page


WORK ORDER #: 1011113A

Work Order Summary

CLIENT:	Mr. Brian Baker Environmental Health & Engineering, Inc. 117 Fourth Avenue Needham, MA 02494	BILL TO:	Accounts Payable Environmental Health & Engineering, Inc. 117 Fourth Avenue Needham, MA 02494
PHONE:	800-825-5343	P.O. #	17131
FAX:	781-247-4305	PROJECT #	17131
DATE RECEIVED:	11/04/2010	CONTACT:	Ausha Scott
DATE COMPLETED:	11/15/2010		

<u>FRACTION #</u>	<u>NAME</u>	<u>TEST</u>
17B	Lab Blank	ATL Applications
18A	LCS	ATL Applications

CERTIFIED BY:



Laboratory Director

DATE: 11/15/10

This report shall not be reproduced, except in full, without the written approval of Air Toxics Ltd.

180 BLUE RAVINE ROAD, SUITE B FOLSOM, CA - 95630
(916) 985-1000 . (800) 985-5955 . FAX (916) 985-1020

**LABORATORY NARRATIVE
Hydrogen Sulfide by Radiello 170
Environmental Health & Engineering, Inc.
Workorder# 101113A**

Sixteen Radiello 170 (H₂S) samples were received on November 04, 2010. The procedure involves adsorption of H₂S by zinc acetate to form zinc sulfide. The sulfide is then recovered by extraction with water and addition of ferric chloride in a strongly acidic solution to produce methylene blue. Methylene blue absorbance is then measured at 665 nm using a spectrophotometer. Results are reported in uG and uG/m³.

Sampling rate of 69 mL/min for H₂S was provided by the manufacturer.

Receiving Notes

There were no receiving discrepancies.

Analytical Notes

Results were calculated based on 25 deg C without temperature correction. The actual exposure time was used to calculate sample concentrations and reporting limits.

An exposure time of 18,700 minutes was used for the QC samples.

Definition of Data Qualifying Flags

Eight qualifiers may have been used on the data analysis sheets and indicate as follows:

- B - Compound present in laboratory blank greater than reporting limit.
- J - Estimated value.
- E - Exceeds instrument calibration range.
- S - Saturated peak.
- Q - Exceeds quality control limits.
- U - Compound analyzed for but not detected above the detection limit.
- M - Reported value may be biased due to apparent matrix interferences.
- N - The identification is based on presumptive evidence.

File extensions may have been used on the data analysis sheets and indicates as follows:

- a-File was requantified
- b-File was quantified by a second column and detector
- r1-File was requantified for the purpose of reissue

Sample Results and Raw Data

AIR TOXICS LTD.

ATL Application # 59 for RAD 170 (Hydrogen Sulfide)

Spectrophotometer

Field Sample I.D.	Lab Sample I.D.	Collection Date	Analysis Date	Dilution Factor	Reporting Limit (ug)	Reporting Limit (ug/m3)	Amount (ug)	Amount (ug/m3)
118581	1011113A-01A	11/1/2010	11/12/2010	1.00	0.80	0.58	ND	ND
118582	1011113A-02A	11/1/2010	11/12/2010	1.00	0.80	0.58	ND	ND
118583	1011113A-03A	11/1/2010	11/12/2010	1.00	0.80	0.58	ND	ND
118584	1011113A-04A	11/1/2010	11/12/2010	1.00	0.80	0.58	1.3	0.94
118585	1011113A-05A	NA	11/12/2010	1.00	0.80	0.58	ND	ND
118586	1011113A-06A	NA	11/12/2010	1.00	0.80	0.58	ND	ND
118597	1011113A-07A	11/1/2010	11/12/2010	1.00	0.80	0.58	1.9	1.4
118598	1011113A-08A	11/1/2010	11/12/2010	1.00	0.80	0.58	1.7	1.3
118599	1011113A-09A	11/1/2010	11/12/2010	1.00	0.80	0.58	1.7	1.3
118600	1011113A-10A	11/1/2010	11/12/2010	1.00	0.80	0.58	ND	ND
118601	1011113A-11A	NA	11/12/2010	1.00	0.80	0.58	ND	ND
118602	1011113A-12A	NA	11/12/2010	1.00	0.80	0.58	ND	ND
118613	1011113A-13A	11/1/2010	11/12/2010	1.00	0.80	0.58	1.3	0.94
118614	1011113A-14A	11/1/2010	11/12/2010	1.00	0.80	0.58	ND	ND
118615	1011113A-15A	11/1/2010	11/12/2010	1.00	0.80	0.58	1.6	1.2
118616	1011113A-16A	11/1/2010	11/12/2010	1.00	0.80	0.58	ND	ND
118616 Lab Duplicate	1011113A-16AA	11/1/2010	11/12/2010	1.00	0.80	0.58	ND	ND
Method Blank	1011113A-17A	NA	11/12/2010	1.00	0.80	0.58	ND	ND
Method Blank	1011113A-17B	NA	11/12/2010	1.00	0.80	0.58	ND	ND
LCS	1011113A-18A	NA	11/12/2010	1.00	0.80	0.58	%Rec 100	

COMMENTS: 1. NA=Not Applicable

2. ND=Not Detected

3. Exposure time of 18700 minutes was assumed for the QC samples.

4. Background subtraction not performed.

Hydrogen Sulfide Radiello Calculation Worksheet

Workorder #: 101113A

Sampling Rate (ng/ppb.min) 0.096 Typically 0.096 for H2S

Sampling T (deg C) 25 Typically 25

Volume (ml) 10.5 Typically 10.5 for H2S

Date of Analysis: 11/12/2010

Corrected Q 0.096 Takes into account temp

(Abs-Y:ln)DXF Slope

Conc(ug/ml)xVol (ml)

conc (ug sulfide) *MW H2S
MW Sulfide

Q includes conversion from Sulfide to H2S

Conc (ug) x 1000
Q x Duration

ppbx mw
24.45

T Corrected, no Blank correction

LabSampleID	Client	Date of Collection	Abs	Duration (min)	DF	Conc (ug/ml) of sulfide	Conc (ug) of sulfide	Conc (ug) of H2S	Conc (ppb) of H2S	Conc (ug/m3) of H2S
01A	118581	11/1/2010	0.069	18695	1.00	0.028575316	0.300040813	0.318864702	0.167	0.238
02A	118582	11/1/2010	0.077	18695	1.00	0.036004943	0.3780519	0.401770029	0.211	0.294
03A	118583	11/1/2010	0.076	18695	1.00	0.035076239	0.368300514	0.391406863	0.205	0.286
04A	118584	11/1/2010	0.163	18695	1.00	0.115873436	1.216671081	1.2930023	0.678	0.945
05A	118585	NA	0.021	18700	1.00	-0.016002448	-0.168025707	-0.178567263	-0.094	-0.130
06A	118586	NA	0.025	18700	1.00	-0.012287635	-0.129020163	-0.137114599	-0.072	-0.100
07A	118597	11/1/2010	0.224	18700	1.00	0.172524344	1.811505616	1.925155423	1.009	1.407
08A	118598	11/1/2010	0.204	18700	1.00	0.153950276	1.6164779	1.717892104	0.900	1.255
09A	118599	11/1/2010	0.207	18700	1.00	0.156736386	1.645732057	1.748981602	0.917	1.278
10A	118600	11/1/2010	0.07	18700	1.00	0.029504019	0.309792199	0.329227868	0.173	0.241
11A	118601	NA	0.022	18700	1.00	-0.015073745	-0.158274321	-0.168204097	-0.088	-0.123
12A	118602	NA	0.021	18700	1.00	-0.016002448	-0.168025707	-0.178567263	-0.094	-0.130
13A	118613	11/1/2010	0.163	18690	1.00	0.115873436	1.216671081	1.2930023	0.678	0.945
14A	118614	11/1/2010	0.095	18690	1.00	0.052721604	0.553576845	0.588307016	0.309	0.430
15A	118615	11/1/2010	0.191	18690	1.00	0.14187132	1.489709884	1.583170947	0.830	1.157
16A	118616	11/1/2010	0.073	18690	1.00	0.032290129	0.339046356	0.360317366	0.189	0.263
16AA	118616 Lab Duplicate	11/1/2010	0.076	18690	1.00	0.035076239	0.368300514	0.391406863	0.205	0.286
					1.00	-0.03550522	-0.372804809	-0.396193748	#DIV/0!	#DIV/0!
					1.00	-0.03550522	-0.372804809	-0.396193748	#DIV/0!	#DIV/0!
					1.00	-0.03550522	-0.372804809	-0.396193748	#DIV/0!	#DIV/0!
					1.00	-0.03550522	-0.372804809	-0.396193748	#DIV/0!	#DIV/0!
					1.00	-0.03550522	-0.372804809	-0.396193748	#DIV/0!	#DIV/0!
17A	Method Blank	NA	0.023	18700	1.00	-0.014145041	-0.148522935	-0.157840931	-0.083	-0.115
17B	Method Blank	NA	0.021	18700	1.00	-0.016002448	-0.168025707	-0.178567263	-0.094	-0.130
18A	LCS	NA	0.181	18700	1.00	0.132590098	1.392196026	1.479539287	0.776	1.081

QC Duration
18700

CCV Spike Amt
0.133

Calibration Date

11/12/2010 Linear Regression

Slope	1.076770029
Y-int	0.038230957
R2	0.998295164

	Slope	Y-Int	R2
0	1.076770029	0.038230957	0.998295164
0.0716	0.098		
0.143	0.187		
0.286	0.354		
0.572	0.684		
1.145	1.256		

QC Results and Raw Data

Work Order: 101113ADate: 11/12/10Method: Rad 170Analyst: M. SkidmoreWavelength: 665nm

Standard ID	Concentration	ABS
	Sulfide (ug/mL)	
Level 1 1993-91-E	0.0716	0.098
Level 2 -D	0.143	0.187
Level 3 -C	0.286	0.354
Level 4 -B	0.572	0.684
Level 5 -A	1.145	1.256
ICV 1993-92	0.286	0.351

$$r = 0.9983$$

$$m = 1.077$$

$$b = 0.03823$$

ICV % Recovery = 102

Fraction	Dilution	ABS	Sample ID	Sample Volume	Comments
01A	1.00	0.069	118581	10.5 mL	
02A		0.077	118582		
03A		0.076	118583		
04A		0.163	118584		
05A		0.021	118585		
06A		0.025	118586		
07A		0.224	118597		
08A		0.204	118598		
09A		0.207	118599		
10A		0.070	118600		
11A		0.022	118601		
12A		0.021	118602		
13A		0.163	118613		
14A		0.095	118614		
15A		0.191	118615		
16A		0.073	118616		
16AA		0.076	118616		
BIK1		0.023	N/A		Lot: 10101
BIK2		0.021			↓
LCS		0.181			Lot: 10101, 0.133 ug/mL
CCV	✓	0.356	↓	5.0 mL	0.286 ug/mL

Procedure:

- 1.) Add 10 mL of H₂O to sample tube, cap and vortex for 1 minute.
- 2.) Add 0.5 mL of Ferric Chloride-Amine solution and cap immediately.
- 3.) Allow color to develop for 30 minutes.
- 4.) Measure absorbance at 665nm.

MJS 11/12/10

Mil. B. B.
Signed

11/12/10
Date

Spectrophotometer Standard Preparation Log

@Air Toxics Ltd.

Log Book #: 1993

Standard ID: 1993-76

Project: Rad 170 Amine Solution

Analyst: M. Skidmore

Preparation Date: 10/18/10

Expiration Date: 11/18/10

Solvent: HPLC H₂O

Solvent Lot #: DB 270

Procedure/Comments:

Sulfuric Acid Solution:

Slowly add 6.25 mL of concentrated sulfuric acid to 2.5 mL of D.I. H₂O, and let the solution cool. (sulfuric acid lot: 01428LS).

Amine Solution:

Dissolve 1.6875g of N,N-dimethyl-p-phenyldiammonium oxalate (located in ER1A; Lot: 63797PJ) in the above mentioned sulfuric acid solution. Dilute this solution to 250 mL with sulfuric acid-water 1:1 v/v. (This is roughly 120 mL H₂O + 120 mL sulfuric acid).

MJS 10/18/10

MJS 10/18/10

Spectrophotometer Standard Preparation Log

@Air Toxics Ltd.

Log Book #: 1993

Standard ID: 1993-77

Project: Ferric Chloride Solution Rad 170

Analyst: M. Skidmore


Preparation Date: 10/18/10

Expiration Date: 10/18/11

Solvent: HPLC H₂O

Solvent Lot #: DB 270

Procedure/Comments: Dissolve 125 g of ferric chloride hexahydrate
(located in ERAC, lot: 73297) in 50 mL of H₂O,



MJS 10/18/10

Miles R. O.

Signed

10/18/10

Date

Fauzi

Reviewed

10/22/10

Date

Spectrophotometer Standard Preparation Log

@Air Toxics Ltd. Log Book #: 1993

Standard ID: 1993-90

Project: Ferric chloride Amine solution

Analyst: M. Skidmore

Preparation Date: 11/12/10

Expiration Date: 11/12/10

Solvent: HPCL H₂O

Solvent Lot #: DB 812

Procedure/Comments: Add 4.0 mL of ferric chloride solution
(1993-77, exp 10/18/11) with 20 mL of amine solution
(1993-76, exp 11/18/10).

MJS 11/12/10

Spectrophotometer Standard Preparation Log

@Air Toxics Ltd. Log Book #: 1993

Standard ID: 1993-91

Project: Rad 170 calibration curve

Analyst: M. Skidmore

Preparation Date: 11/12/10

Expiration Date: 11/12/10

Solvent: HPLC H₂O

Solvent Lot #: DB812

Procedure/Comments: _____

Solution A: 2 mL of Code Rad 171 (1476-2077, exp 6/16/11) (located in ER1B) with 98 mL of D.I. H₂O = 1.145 µg/mL

Solution B: 2.5 mL of Solution A with 2.5 mL of D.I. H₂O = 0.572 µg/mL

Solution C: 1.25 mL of Solution A with 3.75 mL of D.I. H₂O = 0.286 µg/mL

Solution D: 0.625 mL of Solution A with 4.375 mL of D.I. H₂O = 0.143 µg/mL

Solution E: 0.375 mL of Solution A with 5.625 mL of D.I. H₂O = 0.0716 µg/mL

Note: Each solution was measured immediately after it was prepared. Solution A is only stable in the flask it was prepared in.

MJS 11/12/10

MJS
11/12/10

Spectrophotometer Standard Preparation Log

@Air Toxics Ltd. Log Book #: 1993

Standard ID: 1993-92
Project: Rad. 170 ICV
Analyst: Fm
Preparation Date: 11/12/10
Expiration Date: 11/12/10

Solvent: HPLC H2O
Solvent Lot #: DB 812

Procedure/Comments: _____

Solution A: 2 mL of Code Rad 171 (1476-2077, exp 6/16/11) (located in ER1B) with
98 mL of D.I. H₂O = 1.145 µg/mL

Solution C: 1.25 mL of Solution A with 3.75 mL of D.I. H₂O = 0.286 µg/mL

Note: Each solution was measured immediately after it was prepared. Solution A is only
stable in the flask it was prepared in.

Fm 11/12/10

Fm
11/12/10

Fauzi
Signed

11/12/10
Date

Mike B...
Reviewed

11/12/10
Date

Spectrophotometer Standard Preparation Log

@Air Toxics Ltd. Log Book #: 1993


Standard ID: 1993-93

Project: Ferric Chloride-Amine Solution Rad 170⁺
Analyst: M. Skidmore H₂S LCSPreparation Date: 11/12/10Expiration Date: 11/12/10Solvent: MPLC H₂OSolvent Lot #: DB 812

Procedure/Comments:

A Rad 170 cartridge (lot: 10101) was placed in a 40 mL VOA vial. 10.0 mL of D.I. H₂O was aliquoted into the vial. 1.0 mL of H₂S gas (1476-1497, 1000 ppm) was injected into the vial, into the H₂O. The solution was allowed to gently shake for 2 hours. Then 0.5 of the ferric-chloride-amine (1993-90) was added to the vial and capped immediately. The solution was allowed to sit for 30 minutes and the absorbance was measured at 665 nm.

MJS 11/12/10

MJS
11/12/10

Shipping/ Receiving Documents

180 Blue Ravine Road, Suite B
Folsom, CA 95630

Phone (916) 985-1000 FAX (916) 985-1020
Hours 8:00 A.M. to 6:00 P.M. Pacific

COMPANY: Environmental Health & Engineering, Inc.
ATTENTION: Mr. Brian Baker
FAX #: 781-247-4305
FROM: Sample Receiving
Workorder #: 1011113A
of pages (Including Cover): 4
11/16/2010

Thank you for selecting Air Toxics Ltd. We have received your samples and have found discrepancies. In order to expedite analysis and reporting, please review the attached information for accuracy. Corrections can be faxed to **Ausha Scott at 916-985-1020.**

ATL will proceed with the analysis as specified on the Chain of Custody and Sample Login page.

In accordance with your company's contract, this account is required to have a PO that is fully executed by both parties which also covers the cost of the workorder before any data can be released. Please ensure that you have given all appropriate information to our Project Manager so that there will be no delay in reporting of the data you are requesting.

Your prompt response is appreciated.

CHAIN OF CUSTODY FORM

DATE: 3 NOV 10

1011113

FROM: Environmental Health and Engineering, Inc.
117 Fourth Avenue
Needham, MA 02494-2725

TO: AIR TOXICS

Please send invoices to ATTN: Accounts Payable
Please send reports to ATTN: Data Coordinator

In all correspondence regarding this matter, please refer to EH&E Project # 17131

The cost of this analysis will be covered by EH&E Purchase Order # 17131

For EH & E Data Coordinator - URGENT DATA ☐

	SAMPLE ID	SAMPLE TYPE	ANALYTICAL METHOD/NUMBER		OTHER:Time/Date/Vol.
01A	118581	PASSIVE AIR	H ₂ S ANALYSIS	10/19/10 - 11/1/10	12D 23H 35M
02A	118582			I	I
03A	118583			I	I
04A	118584			I	I
05A	118585				60
06A	118586				I
07A	118597			10/19/10 - 11/1/10	12D 23H 40M
08A	118598			I	I
09A	118599			I	I
10A	118600			I	I
11A	118601				60
12A	118602				I
13A	118613			10/19/10 - 11/1/10	12D 23H 30M
14A	118614			I	I
15A	118615			I	I
16A	118616			I	I

Special instructions:

- ☐ Standard turn around time ☐ Rush by _____ date/time ☐ Other _____
☐ Fax results 781-247-4305
☐ RETURN SAMPLES ☒ Electronic transfer - datacoordinator@ehinc.com
☒ Additional report recipient Hakere@ehinc.com

Each signatory please return one copy of this form to the above address

Relinquished by: [Signature] of Environmental Health & Engineering, Inc. Date: 11/3/10

Received by: Frank Whittaker of (company name) ATC Date: 11/4/10 900

Relinquished by: _____ of (company name) Fedex Date: _____

Received by: _____ of (company name) CUSTODY SEAL INTACT? Date: _____

Relinquished by: _____ of (company name) Y N NON-TEMP 3.5C Date: _____

Received by: _____ of (company name) Date: _____

Received by: _____ of Environmental Health & Engineering, Inc. Date: _____

SAMPLE RECEIPT SUMMARY

WORKORDER 1011113A

Client

Mr. Brian Baker
Environmental Health &
Engineering, Inc.
117 Fourth Avenue
Needham, MA 02494

Phone

800-825-5343

Fax

781-247-4305

Date Promised: 11/17/10 11:59 pm

Date Completed: 11/15/10

Date Received: 11/4/10

PO#: 17131

Project#: 17131

Total \$: \$ 1,360.00

Sales Rep: TL

Logged By: MW

<u>Fraction</u>	<u>Sample #</u>	<u>Analysis</u>	<u>Collected</u>	<u>Amount\$</u>
01A	118581	ATL Applications	11/1/2010	\$80.00
02A	118582	ATL Applications	11/1/2010	\$80.00
03A	118583	ATL Applications	11/1/2010	\$80.00
04A	118584	ATL Applications	11/1/2010	\$80.00
05A	118585	ATL Applications	NA	\$80.00
06A	118586	ATL Applications	NA	\$80.00
07A	118597	ATL Applications	11/1/2010	\$80.00
08A	118598	ATL Applications	11/1/2010	\$80.00
09A	118599	ATL Applications	11/1/2010	\$80.00
10A	118600	ATL Applications	11/1/2010	\$80.00
11A	118601	ATL Applications	NA	\$80.00
12A	118602	ATL Applications	NA	\$80.00
13A	118613	ATL Applications	11/1/2010	\$80.00
14A	118614	ATL Applications	11/1/2010	\$80.00
15A	118615	ATL Applications	11/1/2010	\$80.00
16A	118616	ATL Applications	11/1/2010	\$80.00
16AA	118616 Lab Duplicate	ATL Applications	11/1/2010	\$0.00
17A	Lab Blank	ATL Applications	NA	\$0.00
17B	Lab Blank	ATL Applications	NA	\$0.00
18A	LCS	ATL Applications	NA	\$0.00

Note: Samples received after 3 P.M. PST are considered to be received on the following work day.
Atlas Project Name/Profile#: CPSC/14482

BILL TO: Accounts Payable
Environmental Health & Engineering, Inc.
117 Fourth Avenue
Needham, MA 02494

Analysis Code: Other GC

TERMS:

Reporting Method: ATL Application #59 H2S-Radiello 170

180 BLUE RAVINE ROAD, SUITE B FOLSOM, CA - 95630
(916) 985-1000 . (800) 985-5955 . FAX (916) 985-1020

SAMPLE RECEIPT SUMMARY Continued

Client

Phone

Date Promised:

Date Completed:

Date Received:

Fax

PO#:

Project#:

Sales Rep:

Total \$: \$ 1,360.00

Logged By: MW

<u>Fraction</u>	<u>Sample #</u>	<u>Analysis</u>	<u>Collected</u>	<u>Amount\$</u>
Misc. Charges eCVP (16) @ \$5.00 each.				\$80.00

Note: Samples received after 3 P.M. PST are considered to be received on the following work day.
 Atlas Project Name/Profile#: CPSC/14482

BILL TO: Accounts Payable
 Environmental Health & Engineering, Inc.
 117 Fourth Avenue
 Needham, MA 02494

Analysis Code: Other GC

TERMS:

Reporting Method: ATL Application #59 H2S-Radiello 170

180 BLUE RAVINE ROAD, SUITE B FOLSOM, CA - 95630
(916) 985-1000 . (800) 985-5955 . FAX (916) 985-1020

@ Air Toxics Ltd	Title: Sample Discrepancy Report			Release Date: 03/03/10
	Form #: F1.3	Revision #: 1	Revision Date: 10/7/08	Page #: 1 of 2

Sample Discrepancy Report

Identification

Initiated By: MW Project ID: 14482 PM: AS Date: 11/6/2010 Discrepancy Type: ☐ 1. ☒ 2. ☐ 3.

Workorder(s) affected: 1011113 Sample(s) affected: 05A, 06A, 11A, 12A, 17A, 18A, 23A, 24A, 29A and 30A

1. Sample Receipt Discrepancies

Narration Not Required:

- 1.1. ☐ Sample container (cartridge/tube/VOA vial) was received broken, however sample was intact.
- 1.2. ☐ No brass cap on canister.
- 1.3. ☐ Date of Collection noted on first sample, but no arrow down to indicate all samples.

Notify Lab for further determination:

- 1.4. ☐ Tedlar bag received with minimal volume.

Initials: _____ Date: _____

Narration Required in Lab Narrative and Sample Confirmation:

- 1.5. ☐ COC was not filled out in ink.
- 1.6. ☐ COC improperly relinquished / received.
- 1.7. ☐ Sample tags / can numbers do not match the COC.
- 1.8. ☐ Sample date ☐ error / ☐ missing on COC but noted on sample tag (check one).
- 1.9. ☐ Custody Seal on the outside of the container was ☐ broken / ☐ improperly placed (check one).
- 1.10. ☐ ID-none on the sample Tag/Blank
- 1.11. ☐ Other (describe below).

Describe the Discrepancy: _____

2. Sample Receipt/Screening Discrepancies requiring PM notification

Document on Cover Page of Sample Receipt Confirmation and in Receiving Notes of Lab Narrative

If Section II. is filled out PM must be notified within 24 hrs of initiation

- | | |
|---|---|
| <ul style="list-style-type: none"> 2.1. <input type="checkbox"/> COC was not received with samples. 2.2. <input type="checkbox"/> Analysis method(s) is <input type="checkbox"/> not specified / <input type="checkbox"/> incorrectly specified (check one) on the COC. 2.3. <input type="checkbox"/> Incorrect sampling media / container for analysis requested. 2.4. <input type="checkbox"/> Number of samples on the COC does not match the number of samples that were received. 2.5. <input type="checkbox"/> Samples were received expired. 2.6. <input checked="" type="checkbox"/> Sampling date (time for sulfur) is not documented for <input checked="" type="checkbox"/> <u>some</u> / <input type="checkbox"/> <u>any</u> samples (check one). 2.7. <input type="checkbox"/> Sample received with amount of H₂O in the Tedlar Bag. 2.8. <input type="checkbox"/> Sample cannot be analyzed. Container was <input type="checkbox"/> received broken / <input type="checkbox"/> leaking / <input type="checkbox"/> flat / <input type="checkbox"/> defective. 2.9. <input type="checkbox"/> Tedlar bag / canister received emitting a strong odor; Sample <input type="checkbox"/> can / <input type="checkbox"/> cannot (check one) be analyzed. 2.10. <input type="checkbox"/> Tedlar Bag for Sulfur analysis has metal fitting. 2.11. <input type="checkbox"/> Environmental Supply Company valves 2.12. <input type="checkbox"/> Sorbent samples-sampling volume was not provided | <ul style="list-style-type: none"> 2.13. <input type="checkbox"/> Flow controller used – canister samples received at ambient or under pressure. 2.14. <input type="checkbox"/> Canister was at ambient pressure at time of pressurization and (check all that apply):
 <input type="checkbox"/> Canister failed leak check on two manifolds,
 <input type="checkbox"/> Canister valve was open,
 <input type="checkbox"/> Brass nut was loose/not present.
 <input type="checkbox"/> Sample can be analyzed
 <input type="checkbox"/> Cannot be analyzed 2.15. <input type="checkbox"/> Canister sample received with a vacuum difference >5.0"Hg between the receipt vac. And the final vac. reported on the COC, indicating loss of vacuum. 2.16. <input type="checkbox"/> Canister sample received at >15"Hg (<u>not</u> identified as a Trip/Field Blank). 2.17. <input type="checkbox"/> Canister Trip Blank received at low vacuum (< 25"Hg). 2.18. <input type="checkbox"/> Sorbent Sample received outside method required temperature of 2°C to 6°C; <input type="checkbox"/> ice / <input type="checkbox"/> blue ice (check one) was present. A temp. Blank <input type="checkbox"/> was / <input type="checkbox"/> was not present (check one). 2.19. <input type="checkbox"/> Other (describe below) |
|---|---|

Initials: _____ Date: _____ Notify Receiving: ☐ Notify PM: ☐

Describe the Discrepancy: No DOCs for samples listed above (NO NEED TONARRATE)

3. Lab Discrepancies requiring Team Leader/PM notification

Document in Analytical Notes of Lab Narrative

If Section III. is filled out PM must be notified within 24 hrs of initiation

- | | |
|--|--|
| 3.1. <input type="checkbox"/> Tedlar Bag found to be leaking at the time of analysis; sample <input type="checkbox"/> can / <input type="checkbox"/> cannot (check one) be analyzed. | 3.6. <input type="checkbox"/> Sample loss due to instrument malfunction / broken glassware. |
| 3.2. <input type="checkbox"/> Tedlar Bag found to be flat/low volume; sample cannot be analyzed. | 3.7. <input type="checkbox"/> Low/high surrogate recoveries noted in QC/sample(s) for extractable samples. |
| 3.3. <input type="checkbox"/> Sulfur samples received with insufficient time to analyze prior to expiration. | 3.8. <input type="checkbox"/> Reporting Limit was raised. |
| 3.4. <input type="checkbox"/> Canister found to be leaking at the time of analysis. | 3.9. <input type="checkbox"/> Post weight > Pre weight in field/lab Blank for PM10/TSP samples. |
| 3.5. <input type="checkbox"/> VOST tube saturated; bag dilution necessary. | 3.10. <input type="checkbox"/> Other (describe below). |

Initials: _____ Date: _____ Notify Receiving: ☐ Notify PM: ☐

Team Lead Initials: _____ Date: _____

Describe the Discrepancy: _____

How Does this Affect Client: _____

Project Manager Use Only

Project Manager Notification

☒ Section 2 Complete

☐ Section 3 Complete

Action:

- ☒ It is not necessary to notify the client. Narrate the discrepancy in Receiving Notes/Analytical Notes of Lab Narrative.

PM Initials: AS Date: 11/8/2010

- ☐ Client notification required. See attached client contact / email, or comments below:

Client Notification:

PM Initials: _____ Person notified: _____ Date: _____

- ☐ Waiting for Client Reply

Samplers that have "0D" are trip blanks. They will not have DOCs. Please do not discrepancy these.

Comments: _____

☐ Notify Lab Name: _____ Date: _____ Notify Receiving: ☒

- ☐ Additional notifications attached.

Additional Comments:

Other Records

Method : ATL Application #59 H2S-Radiello 170

CAS Number	Compound	Rpt. Limit (ug)
7783-06-4	Hydrogen Sulfide	1.2

@ Air Toxics Ltd	Title: Data Review Checklist		Release Date: 07/28/10	
	Form #: F1.27	Revision #: 2	Revision Date: 07/27/10	Page #: 1 of 2

DATA REVIEW CHECKLIST

Work Order #:

101113A

A ₁	A ₂	W	T	R	Q	
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Analysis/Reporting vs. Project Profile/SOP requirements checked (i.e. 100% Dups, J-Flag to MDL, etc)
		<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	The final report has the correct reporting list, special units, and header info.
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Non-Standard sublist printed/verified, LOQ and LOD verified
		<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Lab Narrative is correct (proper method & description/Receiving & Analytical notes correct)
		<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Sample Discrepancy Report (SDR) is completed
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Corrective Action issued - #
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Unusual circumstances have been documented in the notes section below
						LUMEN validation report present and initialed
						CIRCLE (YES / NO)
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Lab Blank, CCV, LCS and DUP met QC criteria
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Hold time is met for all samples
		<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Appropriate data qualifier flags are applied
		<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Manual integrations for samples and QC are properly documented
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Samples analyzed within the project or method specific clock
		<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Retention times have been verified
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Appropriate ICAL(s) included, %RSD Recalculation
		<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	At least one result per sample is verified against the target quant sheets/raw data
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Dilution factor correctly calculated (sample load volume, syringe and bag dilutions, can pressurization(s))
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Correct amount of sample analyzed (i.e. sample not over-diluted)
		<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Spectra verified - documentation of spectral defense included (Section 5A of eCVP pkg)
		<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	TICs resemble reference spectra
		<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	TICs between duplicate samples are consistent
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Checked samples for trends (i.e. Influent vs. Effluent, Field Dups, Field/Trip Blank, etc.)
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Data for multiple analyses of sample(s) has been evaluated for comparability of results
		<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Special units for all samples in the final report are correctly calculated
		<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Manually entered results checked (i.e. TPH/NMOC)
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Chain of Custody verified for any special comments (i.e. different compounds/RLs, action levels)
		<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Chain of Custody scanned correctly
		<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Verify sample id's vs. chain of custody
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Date MDL(s) performed per instrument(s) 10/25/10
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Samples pressurized w/ appropriate gas (N ₂ or He) <input type="checkbox"/> Other (i.e. Tedlar bag, cartridge, sorbent)
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Final pressure consistent with canister size (6L vs. 1L)
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Verify receipt pressures
<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Verify canister ID #'s
		<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Final invoice amount correct (adjusted for TAT, Penalties, Re-issue Charges etc.)
		<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Final PDF report reviewed for correctness

Notes: (to include: noting samples with QA/QC problems, Blanks with positive hits, narratives, etc.)

A/R: 16A - Dup

T/Q:

A ₁ /A ₂	W/T	R*	Q
(Analytical Review/Date)	(Write-up/Tech Review/Date)	(Report Review/Date)	(QA Review/Date)
A ₁ : <i>[Signature]</i> 11/12/10	W: <i>[Signature]</i> 11/12/10	R:	

A₂:

T:

Note (1): Please check all the appropriate boxes. Indicate "NA" for any statement that does not apply.

Note (2): Report reviewer and write-up reviewer must be separate individuals for DoD & Client Specific projects.

* Report Review is completed for DoD & Client Specific projects only.